

Package ‘AmpliconDuo’

January 20, 2025

Type Package

Title Statistical Analysis of Amplicon Data of the Same Sample to Identify Artefacts

Version 1.1.1

Date 2020-05-22

Author Anja Lange [aut, cre],
Daniel Hoffmann [aut]

Maintainer Anja Lange <anja.lange@uni-due.de>

Depends R (>= 2.10), stats, ggplot2, xtable

Description Increasingly powerful techniques for high-throughput sequencing open the possibility to comprehensively characterize microbial communities, including rare species. However, a still unresolved issue are the substantial error rates in the experimental process generating these sequences. To overcome these limitations we propose an approach, where each sample is split and the same amplification and sequencing protocol is applied to both halves. This procedure should allow to detect likely PCR and sequencing artifacts, and true rare species by comparison of the results of both parts. The AmpliconDuo package, whereas amplicon duo from here on refers to the two amplicon data sets of a split sample, is intended to help interpret the obtained read frequency distribution across split samples, and to filter the false positive reads.

License GPL (>= 3)

NeedsCompilation no

Repository CRAN

Date/Publication 2020-05-25 22:20:02 UTC

Contents

AmpliconDuo-package	2
accepted.amplicons	4
ampliconduo	5
ampliconfreqs	7
ampliconfreqs.long	8
amplicons	9

discordance.delta	9
filter.ampliconduo	11
filter.ampliconduo.set	12
plotAmpliconduo	13
plotAmpliconduo.set	15
plotORDensity	17
site.f	19

Index	20
--------------	-----------

AmpliconDuo-package	<i>Statistical Analysis Of Amplicon Data Of The Same Sample To Identify Spurious Amplicons</i>
---------------------	--

Description

Increasingly powerful techniques for high-throughput sequencing open the possibility to comprehensively characterize microbial communities, including rare species. However, a still unresolved issue are the substantial error rates in the experimental process generating these sequences. To overcome these limitations we propose an approach, where each sample is split and the same amplification and sequencing protocol is applied to both halves. This procedure should allow to detect likely PCR and sequencing artifacts, and true rare species by comparison of the results of both parts.

The **AmpliconDuo** package, whereas ampliconduo from here on refers to the two amplicon data sets of a split sample, is intended to help interpret the obtained amplicon frequency distribution across split samples, and to filter the false positive amplicons.

Details

Package: AmpliconDuo
 Type: Package
 Version: 1.1.1
 Date: 2020-05-22
 License: GPL-2

The core of this package is the `ampliconduo` function, that generates for each pair of a split samples an ampliconduo data frame, while statistically analysing the data by Fisher's exact test. Ampliconduo data frames, or lists of these, are the input required for all other functions of this package.

- `plotAmpliconduo` plots for an ampliconduo the amplicon frequencies (number of reads per amplicon) of sample A vs. amplicon frequencies of sample B, highlighting amplicons displaying a significant deviation between both samples.
- `plotAmpliconduo.set` does the same as `plotAmpliconduo` but accepts a list of ampliconduo data frames and arranges the plots in a 2-dimensional array.
- `plotORDensity` generates a histogram plot of the amplicon frequency odds ratio density for an ampliconduo data frame. For multiple data frames organizes the plots in a 2-dimensional array.

- `discordance.delta` calculates delta (Δ) and delta prime (Δ'), the fraction of amplicon frequencies and amplicons, respectively, with a false discovery rate below a certain threshold θ as a measure of discordance between two amplicon data sets A and B.
- `filter.ampliconduo` applies filter criteria to an ampliconduo data frame deciding which amplicons are going to be rejected.
- `filter.ampliconduo.set` same as `filter.ampliconduo` for a list of ampliconduo data frames.
- `accepted.amplicons` returns the indices of those amplicons that have passed the filter criteria.

Author(s)

Anja Lange (<anja.lange@uni-due.de>) and Daniel Hoffmann (<daniel.hoffmann@uni-due.de>)
Maintainer: Anja Lange (<anja.lange@uni-due.de>)

References

Lange A, Jost S, Heider D, Bock C, Budeus B, et al. (2015) AmpliconDuo: A Split-Sample Filtering Protocol for High-Throughput Amplicon Sequencing of Microbial Communities. PLOS ONE 10(11): e0141590

Examples

```
## load test amplicon frequency data ampliconfreqs and vector with sample names site.f
data(ampliconfreqs)
data(site.f)

## generating ampliconduo data frames
## depending on the size of the data sets, may take some time
ampliconduo.set <- ampliconduo(ampliconfreqs[,1:4], sample.names = site.f[1:2])

## plot amplicon read numbers of sample A vs. amplicon read numbers of sample B,
## indicating amplicons with significant deviations in their occurrence across samples
plotAmpliconduo.set(ampliconduo.set, nrow = 3)

## calculate discordance between the two data sets of an ampliconduo
discordance <- discordance.delta(ampliconduo.set)

## plot the odds ratio density of ampliconduo data
plotORDensity(ampliconduo.set)

## apply filter criteria to remove/mark spurious amplicons
ampliconduo.set.f <- filter.ampliconduo.set(ampliconduo.set, min.freq = 1, q = 0.05)

## return indices of accepted amplicons, indices correspond to indices of the ampliconfreqs data,
## that were used as input for the ampliconduo function
accep.reads <- accepted.amplicons(ampliconduo.set.f)
```

accepted.amplicons *List Of Amplicon Indices That Passed Filter Criteria*

Description

Returns the indices of those amplicons in an ampliconduo data frame, that passed the applied filter criteria (ampliconduo data frames are filtered using the [filter.ampliconduo](#) or [filter.ampliconduo.set](#) function).

Usage

```
accepted.amplicons(x)
```

Arguments

x An ampliconduo data frame or a list of ampliconduo data frames.

Details

Calling this function on an ampliconduo data frame, or a list of the latter, returns the indices of amplicons that passed the applied filter criteria. For each ampliconduo data frame an integer vector is created, and if a list of ampliconduo data frames is supplied with x, these are pooled in a list. The returned indices correspond to the data originally used to generate the ampliconduo data frames (parameter A and B in the [ampliconduo](#) function call).

Value

If x is an ampliconduo data frame, an integer vector is returned.

In case x is a list of ampliconduo data frames, a list of integer vectors is returned, one for each data frame.

Author(s)

Anja Lange & Daniel Hoffmann

See Also

[filter.ampliconduo](#) and [filter.ampliconduo.set](#)

Examples

```
## load example data
data(amplicons)

## apply filter criteria
ampliconduos.f <- filter.ampliconduo.set(amplicons, q = 0.05)

## return a list with accepted amplicons
good.reads <- accepted.amplicons(ampliconduos.f)
```

ampliconduo *Apply Fisher's Exact Tests To Two Amplicon Frequency Sets Of The Same Sample*

Description

Implements Fisher's exact test to detect amplicons with significant deviating read numbers between two amplicon sets of the same sample. The p-values of the Fisher's exact test are corrected for multiple testing by computation of the false discovery rates q . This function is intended to help identifying reads that may be the results of experimental artefacts. (The calculation can take some time depending on the size of the data sets and the computing power.)

Usage

```
ampliconduo(A, B = NULL, sample.names = NULL, correction = "fdr", ...)
```

Arguments

A	A list or a data frame containing amplicon occurrences / number of reads per amplicon (integer values).
B	Optional. A list or a data frame containing amplicon occurrences.
sample.names	Optional. A vector or list of characters with names for the amplicon pairs.
correction	Optional. Specifies the correction method for the p-values from Fisher's exact test. Accepts one of the following characters: "holm", "hochberg", "hommel", "bonferroni", "BH", "BY", "fdr" and "none". Default is "fdr". For more details see p.adjust .
...	Arguments passed to the internally called fisher.test function.

Details

If only A is specified, it is assumed that the list elements 1 & 2, 3 & 4 etc. of A are amplicon data of the same sample. In case A and B are specified, the i th frequency set of A and B are combined. For each amplicon data pair, frequencies at the corresponding positions in the lists are assumed to belong to the same amplicon. It is required, that two frequency sets that belong to the same sample, an ampliconduo, have the same length. The ampliconduo function iterates over all amplicon pairs and performs the following tasks:

- amplicons with frequency zero in both samples are removed. Position information is retained.
- For each amplicon Fisher's exact test using the method [fisher.test](#) is performed. The p-value, odds ratio and confidence interval are returned. Via the ..., arguments `conf.level`, `or` and `alternative` can be passed to the [fisher.test](#) function call. Default values are `conf.level = 0.95`, `or = 1` and `alternative = "two.sided"`.
- The p-values are corrected using the [p.adjust](#) function. By default the method by Benjamini & Hochberg (1995) is used. Setting the correction argument to any of the following characters "holm", "hochberg", "hommel", "bonferroni", "BH", "BY", "fdr", "none", the adjustment method for the p-values can be changed. See function [p.adjust](#).

The **AmpliconDuo** package implements further methods to visualize and filter the returned ampliconduo data frames.

Value

A list of data frames, one for each amplicon pair, that will be called *ampliconduo data frame* in the following. List entries are named according to the specified `sample.names` or numbered.

Each ampliconduo data frame has 9 columns

- `freqA`: frequencies of amplicon set A
- `freqB`: frequencies of amplicon set B (taken from argument B if specified)
- `p`: p-values calculated with Fisher's exact test
- `OR`: odds ratio calculated with Fisher's exact test
- `CI.low`: lower confidence limit for OR
- `CI.up`: upper confidence limit for OR
- `rejected`: logical, indicating whether the amplicon was rejected
- `sample`: sample name taken from `sample.name` if specified, same for all rows in a given data frame

Author(s)

Anja Lange and Daniel Hoffmann

References

Y Benjamini and Y Hochberg. Controlling the False Discovery Rate: A Practical and Powerful Approach to Multiple Testing. *Journal of the Royal Statistical Society. Series B (Methodological)*, **57(1)**:289-300, 1995.

See Also

[fisher.test](#), used to calculate the p-value, odds ratio and confidence interval;

[p.adjust](#), called to correct the p-values;

methods to visualize or further manipulate the ampliconduo data frames:

[plotAmpliconduo.set](#), [plotAmpliconduo](#), [discordance.delta](#),

Examples

```
## loads read numbers from example amplicon data sets
data(ampliconfreqs)
data(site.f)

## generate ampliconduo data frames
ampliconduos.a <- ampliconduo(A = ampliconfreqs[,1:4], sample.names = site.f[1:2])
ampliconduos.b <- ampliconduo(A = ampliconfreqs[c(1,3)],
B = ampliconfreqs[c(2,4)], sample.names = site.f[1:2],
conf.level = 0.9)
```

```
## frequency plot  
plotAmpliconduo.set(ampliconduos.a)
```

ampliconfreqs	<i>SSU Amplicon Frequencies Of Different Sampling Sites (Short Version)</i>
---------------	---

Description

A data frame with 16 amplicon data sets taken from 8 different sampling sites. A reduced version (frequencies of 2500 amplicons) of the [ampliconfreqs.long](#) data provided with this package. Samples from each sampling site were sequenced twice, corresponding to sets A and B (e.g. column FU25.A and FU25.B). Names of the sampling sites are specified in the data [site.f](#).

Usage

```
data(ampliconfreqs)
```

Format

A data frame with 80903 observations on the following 16 variables.

FU25.A a numeric vector
FU25.B a numeric vector
FU28.A a numeric vector
FU28.B a numeric vector
FU31.A a numeric vector
FU31.B a numeric vector
FU31.C a numeric vector
FU31.D a numeric vector
FU34.A a numeric vector
FU34.B a numeric vector
FU37.A a numeric vector
FU37.B a numeric vector
UniPond.A a numeric vector
UniPond.B a numeric vector
BogSoil.A a numeric vector
BogSoil.B a numeric vector

References

Boenigk J, Heider D, Jost S, Lange A, Budeus B, Schilling E, Strittmatter A, Hoffmann D: A high-throughput amplicon sequencing and analysis protocol for comparative analyses of microbial communities (*submitted*)

Examples

```
data(ampliconfreqs)
data(site.f)
ampliconduo.a <- ampliconduo(ampliconfreqs[,1:4], sample.names = site.f[1:2])
```

ampliconfreqs.long *SSU Amplicon Frequencies Of Different Sampling Sites*

Description

A data frame with 16 amplicon data sets taken from 8 different sampling sites. Samples from each sampling site were sequenced twice, corresponding to sets A and B (e.g. column FU25.A and FU25.B). Names of the sampling sites are specified in the data [site.f](#).

Usage

```
data(ampliconfreqs.long)
```

Format

A data frame with 80903 observations on the following 16 variables.

FU25.A a numeric vector
FU25.B a numeric vector
FU28.A a numeric vector
FU28.B a numeric vector
FU31.A a numeric vector
FU31.B a numeric vector
FU31.C a numeric vector
FU31.D a numeric vector
FU34.A a numeric vector
FU34.B a numeric vector
FU37.A a numeric vector
FU37.B a numeric vector
UniPond.A a numeric vector
UniPond.B a numeric vector
BogSoil.A a numeric vector
BogSoil.B a numeric vector

References

Boenigk J, Heider D, Jost S, Lange A, Budeus B, Schilling E, Strittmatter A, Hoffmann D: A high-throughput amplicon sequencing and analysis protocol for comparative analyses of microbial communities (*submitted*)

amplicons

Example List of Ampliconduo Data Frames

Description

A list of ampliconduo data frames that was generated calling the `ampliconduo` function with its default parameters and the `ampliconfreqs.long` data provided with this package as input.

Usage

```
data(amplicons)
```

References

Boenigk J, Heider D, Jost S, Lange A, Budeus B, Schilling E, Strittmatter A, Hoffmann D: A high-throughput amplicon sequencing and analysis protocol for comparative analyses of microbial communities (*submitted*)

Examples

```
data(amplicons)
plotAmpliconduo.set(amplicons, nrow = 3)
```

discordance.delta

Measure Of Discordance Between Two Amplicon Data Sets

Description

Calculates delta (Δ) and delta prime (Δ'), the fraction of amplicon frequencies and amplicons, respectively, with a false discovery rate below a certain threshold θ as a measure of discordance between two amplicon data sets A and B.

Usage

```
discordance.delta(x, names = NULL, theta = 0.05, corrected = TRUE,
printToTex = FALSE, directory = NULL, file.name = NULL)
```

Arguments

x	A list of amplicon duo data frames as returned by the <code>ampliconduo</code> function.
names	Optional. Vector or list of characters specifying the sample/amplicon pair names. By default names are taken from the element names of x.
theta	Optional. Numeric, threshold for the false discovery rate. Default value is 0.05.
corrected	Optional. Logical, indicates whether the p-value from Fisher's exact test (FALSE) or the adjusted p-value (TRUE), here called q, is used for calculation of Δ and Δ' .

printToTex	Optional. Logical, if TRUE writes a file with the resulting table in Tex format. Default is FALSE.
directory	Optional. If printToTex is TRUE, the output file will be saved in the specified directory.
file.name	Optional. If printToTex is TRUE, the output file will have the name file.name.

Details

Calculates Δ and Δ' , the fraction of frequencies of amplicons and amplicons, respectively, with false discovery rate below a certain threshold θ as a measure of discordance between two amplicon data sets A and B with occurrence $r_{(iAS)}$ of amplicon i in amplicon set A of sample S (the ampliconduo data frame). Δ and Δ' are defined as follows:

$$\Delta_{S\theta} = \frac{\sum_{i=1}^{n_S} (r_{iAS} + r_{iBS}) \delta(q_{iS} < \theta)}{\sum_{i=1}^{n_S} (r_{iAS} + r_{iBS})}$$

$$\Delta'_{S\theta} = \frac{\sum_{i=1}^{n_S} \delta(q_{iS} < \theta)}{n_S},$$

$$\text{with } \delta(q_{iS} < \theta) = \begin{cases} 1 & \text{for } q_{iS} < \theta \\ 0 & \text{for } q_{iS} \geq \theta \end{cases}$$

for number n_S of amplicons detected in sample/ampliconduo S. $\Delta_{S\theta}$ and $\Delta'_{S\theta}$ are located between 0 (no discordance, i.e. no statistically significant deviations between experimental branches) and 1 (complete discordance).

x is the return value of a ampliconduo call.

Value

Data frame with three columns. The first column contains the sample/amplicon pair names. Second and third column harbor the corresponding Δ and Δ' values, respectively.

Author(s)

Anja Lange & Daniel Hoffmann

See Also

[ampliconduo](#), generates the expected data format for x
[xtable](#), used to convert the returned data frame into a Latex table.

Examples

```
## load example ampliconduo data frame
data(amplicons)

## calculate the discordance between amplicon data sets of an ampliconduo data frame
dd.a <- discordance.delta(amplicons)
dd.b <- discordance.delta(amplicons, theta = 0.1)
```

filter.ampliconduo *AmpliconDuo Amplicon Filter*

Description

Marks or removes amplicons from an ampliconduo data frame according to the specified filter criteria.

Usage

```
filter.ampliconduo(x, min.freq = 1, OR = NULL, q = NULL, p = NULL, remove = FALSE)
```

Arguments

x	Data frame, an ampliconduo data frame returned by the ampliconduo function.
min.freq	Optional. Integer, minimum frequency/read count for an amplicon in each of the two amplicon sets to be retained. Default value is 1.
OR	Optional. Numeric, minimum odds ratio for an amplicon to be retained. If no value is specified the odds ratio is excluded from the filter criteria.
q	Optional. Numeric, minimum value for q, the adjusted p-value for an amplicon to pass the filter. If no value is specified, q is excluded from the filter criteria.
p	Optional. Numeric, minimum p-value for an amplicon to pass the filter. If no value for p is specified p is excluded from the filter criteria.
remove	Optional. Logical, decides whether amplicons that fail the filter criteria should be removed (TRUE), or retained (FALSE).

Details

Takes the ampliconduo data frame *x* and applies each filter criterion that is selected to each amplicon. If an amplicon *i* does not pass each of the applied criteria, the logical value in column `rejected` in row *i* is set to TRUE. In case the parameter `remove` was set to TRUE, all amplicons with `rejected = TRUE` are removed. The position information with respect to the data used as input for the [ampliconduo](#) call is retained.

Value

Data frame corresponding to the input *x*, but with the adjustments in the `rejected` column according to the specified filter criteria, or removed rows (`removed = TRUE`).

Author(s)

Anja Lange & Daniel Hoffmann

See Also

[ampliconduo](#), generates the input data *x* for this method.

[accepted.amplicons](#), returns the indices of amplicons that have passed the filter criteria.

Examples

```
##load example data
data(amplicons)

## extract the first ampliconduo data frame
ampliconduo1 <- amplicons[[1]]

## apply filter criteria
ampliconduo1.f <- filter.ampliconduo(ampliconduo1)
ampliconduo1.f <- filter.ampliconduo(ampliconduo1, min.freq = 2, remove = TRUE)

## to return a list with the indices (corresponding to the indices of the data
## the ampliconduo function was called on) of all amplicons that passed the filter criteria
good.reads <- accepted.amplicons(ampliconduo1.f)
```

```
filter.ampliconduo.set
```

AmpliconDuo List Amplicon Filter

Description

Marks or removes amplicons from each ampliconduo data frame in a list according to the specified filter criteria.

Usage

```
filter.ampliconduo.set(x, min.freq = 1, OR = NULL, q = NULL, p = NULL, remove = FALSE)
```

Arguments

x	List of ampliconduo data frames, return value of an ampliconduo call.
min.freq	Optional. Integer, minimum frequency/read count for a given amplicon in each of the two amplicon sets of an ampliconduo to be retained. Default value is 1.
OR	Optional. Numeric, minimum odds ratio for an amplicon to be retained. If no value is specified the odds ratio is excluded from the filter criteria.
q	Optional. Numeric, minimum value for q, the adjusted p-value for an amplicon to pass the filter. If no value is specified, q is excluded from the filter criteria.
p	Optional. Numeric, minimum p-value for an amplicon to pass the filter. If no value for p is specified p is excluded from the filter criteria.
remove	Optional. Logical, decides whether amplicons that fail the filter criteria should be removed (TRUE), or retained (FALSE).

Details

For every ampliconduo data frame in argument `x`, applies each filter criterion that was specified to each amplicon. If an amplicon i fails any of the applied criteria, the logical value in column `rejected` in row i is set to `TRUE`. In case the parameter `remove` was set to `TRUE`, all amplicons with `rejected = TRUE` are removed. The position information in respect to the data used as input for the `ampliconduo` call are kept. This method uses the function `filter.ampliconduo`.

Value

List of ampliconduo data frames. Same as input parameter `x` but with the adjustments in the `rejected` column according to the specified filter criteria, or removed rows (`removed = TRUE`)

Author(s)

Anja Lange & Daniel Hoffmann

See Also

`filter.ampliconduo`, performs filtering on single ampliconduo data.frames, is called by this method.

`ampliconduo`, generates the input data `x` for this method.

`accepted.amplicons`, returns the indices of amplicons that have passed the filter criteria.

Examples

```
## load example data
data(amplicons)

## apply filter criteria
ampliconduos.f <- filter.ampliconduo.set(amplicons)
ampliconduos.f <- filter.ampliconduo.set(amplicons, min.freq = 3, remove = TRUE)

## to return a list with the indices (corresponding to the indices of the data
## the ampliconduo function was called on) of all amplicons that passed the filter criteria
good.reads <- accepted.amplicons(ampliconduos.f)
```

plotAmpliconduo

Frequency-Frequency Plot Focusing On Significant Deviations In Amplicon Occurrences

Description

Applied to an ampliconduo data frame, one element of the return value of the `ampliconduo` function. Generates a plot of `freqB` over `freqA` (the read numbers of the same amplicon in both halves A and B of a split sample). For amplicons that have significantly deviating read numbers, i.e. with a p-value or adjusted p-value below a certain threshold, points are colored differently (default: red).

Usage

```
plotAmpliconduo(x, color.threshold = 0.05, xlab = "Abundance (PCR A)",
  ylab = "Abundance (PCR B)", main = NULL, log = "xy", corrected = TRUE, asp = 1,
  legend.position = NULL, save = FALSE, path = NULL, file.name = NULL,
  format = "jpeg", h.start = 0, ...)
```

Arguments

x	Ampliconduo data frame, an element of the returned list of the ampliconduo function.
color.threshold	Optional. Numeric value specifying at which p-value or adjusted p-value points in the plot are drawn in complementary color. Default value is 0.05.
xlab	Optional. Character indicating the x-axis label. Default is "Abundance (PCR A)".
ylab	Optional. Character indicating the y-axis label. Default is "Abundance (PCR B)".
main	Optional. Character specifying the overall title of the plot. If no value is passed, takes the sample name from the x\$sample column.
log	Optional. Character specifying the variables to transform to log ("","x", "y", or "xy"). Default is "xy".
corrected	Optional. Logical to decide whether the p-value (FALSE), or corrected p-value (TRUE) is used for coloring. Default value is TRUE.
asp	Optional. Numeric value, the y/x aspect ratio. Default is 1.
legend.position	Optional. Numeric vector of length two. Defines the position of the legend. By default tries to find a position that fits best the arrangement of the plots.
save	Optional. Logical value indicating if the plot should be saved to file. Default value is FALSE.
path	Optional. Character, in case the argument save was set to TRUE, specifies the directory (no backslash or slash at the end) for saving. By default (parameter save is set to TRUE) the plot is saved in the working directory.
file.name	Optional. If argument save was set to TRUE specifies the file name for the plot. The default name is "ampliconduo_<date>.<file extension>"
format	Optional. Character specifying the format of the saved file. One of "eps", "ps", "tex", "pdf", "jpeg", "tiff", "png", "bmp", "svg" and "wmf" (windows only). Default format is "jpeg".
h.start	Optional. Numeric value between 0 and 360 defines the color of the plotted points. Default value is 0 (blue-green, red).
...	Optional. Allows to pass other aesthetics.

Author(s)

Anja Lange & Daniel Hoffmann

References

<https://ggplot2.tidyverse.org/reference/qplot.html>

See Also

ggplot2 package

[qplot](#) internally used to create the plot.

[plotAmpliconduo.set](#), generates a very similar plot for a list of ampliconduo data frames.

[ampliconduo](#), generates the input data.

Examples

```
## load example data
data(amplicons)

## extract the second ampliconduo data frame
ampliconduo2 <- amplicons[[2]]

## plot the amplicon frequencies of the ampliconduo data frame
plotAmpliconduo(ampliconduo2, main = "ampliconduo_2")
plotAmpliconduo(ampliconduo2, main = "ampliconduo_2", h.start = 50, log = "")
plotAmpliconduo(ampliconduo2, h.start = 50, log = "", asp = 2, corrected = FALSE)
```

plotAmpliconduo.set *Frequeny-Frequency Plot Arrangement Focusing On Significant Deviations In Amplicon Occurences*

Description

Called on the return value of the [ampliconduo](#) function, a list of ampliconduo data frames. Generates for each ampliconduo data frame a plot with freqB over freqA and arranges them in a 2-dimensional array, whereas plots in the same row and column share the same scale. Points with a p-value or adjusted p-value below a certain treshold are colored differently (default: red) indicating significant deviations of amplicon occurences between the two samples in an ampliconduo data frame.

Usage

```
plotAmpliconduo.set(x, color.treshold = 0.05, xlab = "Abundance (PCR A)",
  ylab = "Abundance (PCR B)", log = "xy", corrected = TRUE, asp = 1, nrow = 1,
  legend.position = NULL, save = FALSE, path = NULL, file.name = NULL,
  format = "jpeg", h.start = 0, ...)
```

Arguments

x	List of ampliconduo data frames, return value of the ampliconduo function.
color.threshold	Optional. Numeric value specifying at which p-value or adjusted p-value points in the plot are drawn in complementary color. Default value is 0.05.
xlab	Optional. Character indicating the x-axis label. Default is “Abundance (PCR A)”.
ylab	Optional. Character indicating the y-axis label. Default is “Abundance (PCR B)”.
log	Optional. Character specifying the variables to transform to log (“”, “x”, “y”, or “xy”). Default is “xy”.
corrected	Optional. Logical to decide whether the p-value (FALSE), or corrected p-value (TRUE) is used for coloring. Default value is TRUE.
asp	Optional. Numeric value, the y/x aspect ratio. Default is 1.
nrow	Optional. Integer value specifying the number of rows used to arrange the plots. Default is 1.
legend.position	Optional. Numeric vector of length two. Defines the position of the legend. By default tries to find a position that fits the arrangement of the plots best.
save	Optional. Logical value indicating if the plot should be saved to file. Default value is FALSE.
path	Optional. Character, in case the argument save was set to TRUE, specifies the directory (no backslash or slash at the end) for saving. By default (parameter save is set to TRUE) the plot is saved in the working directory.
file.name	Optional. If argument save was set to TRUE specifies the file name for the plot. The default name is “ampliconduo_<date>.<file extension>”
format	Optional. Character specifying the format of the saved file. One of “eps”, “ps”, “tex”, “pdf”, “jpeg”, “tiff”, “png”, “bmp”, “svg” and “wmf” (windows only). Default format is “jpeg”.
h.start	Optional. Numeric value between 0 and 360, defines the color of the plotted points. Default value is 0 (blue-green, red).
...	Optional. Allows to pass other aesthetics.

Details

Generates an arrangement of plots from the return value of the [ampliconduo](#) function, that nicely visualizes those amplicons with a significant deviations in read numbers between the two amplicon data sets. The data in x are transformed and passed to the [qplot](#) function. The 2-dimensional arrangement of the different plots is achieved using [facet_wrap](#). Important aesthetic parameters like color, aspect ratio, legend position ... are easily customized. Optionally, the plot can be saved in a variety of formats.

Author(s)

Anja Lange & Daniel Hoffmann

References

<https://ggplot2.tidyverse.org/reference/qplot.html>

See Also

ggplot2 package

`qplot`, used by `plotAmpliconduo.set` to create the plots.

`facet_wrap`, called for 2-dimensional arrangement of the plots.

`plotAmpliconduo`, generates a very similar plot for a single ampliconduo data frame.

`ampliconduo`, generates the input data, an ampliconduo data frame.

Examples

```
## loads example data of ampliconduo data frames
data(amplicons)

## plot amplicon frequencies of multiple ampliconduo data frames
plotAmpliconduo.set(amplicons[1:4], nrow = 3, h.start = 100)

plotAmpliconduo.set(amplicons[1:4], nrow = 1, corrected = FALSE, color.treshold = 0.1)
```

plotORDensity

Probability Density Plot Of Odds Ratios

Description

Plots for an ampliconduo data frame probability densities of the odds ratios of amplicon occurrences in the two amplicon data sets. The function allows to shift the two extrema (odds ratios OR = 0 and OR = infinity) to the edges of the plot. Plots of multiple ampliconduo data frames are arranged in a 2-dimensional array with shared scales.

Usage

```
plotORDensity(x, log = "x", ncol = 2, adjust.zeroinf = TRUE, zero.pos = 0.005,
  inf.pos = 200, binwidth = 0.15, color = "black", xlab = "odds ratio",
  save = FALSE, path = NULL, file.name = NULL, format = "jpeg", ...)
```

Arguments

x	List or a single ampliconduo data frame, return value of the <code>ampliconduo</code> function.
log	Optional. Character specifying the variables to transform to log (“”, “x”, “y”, or “xy”). Default is “x”.
ncol	Optional. Integer value specifying the number of columns used to arrange the plots. Default is 2.

adjust.zeroinf	Optional. Logical, specifies whether the density bar for 0 and inf should be shifted. Default value is TRUE.
zero.pos	Optional. Numeric, in case adjust.zeroinf is set to TRUE, determines the position for the odds ratio = 0 bar. Default value is 0.05.
inf.pos	Optional. Numeric, in case adjust.zeroinf is set to TRUE, gives the position for the odds ratio = inf bar. Default value is 200.
binwidth	Optional. Numeric, bin width to use, default is 0.15.
color	Optional. Character, name of the color used to draw the density bars. Default is "black".
xlab	Optional. Character, label for the x-axis. Default is "odds ratio".
save	Optional. Logical, TRUE to save the plot, else FALSE. Default is FALSE.
path	Optional. Character, in case the argument save was set to TRUE, specifies the directory (no backslash or slash at the end) for saving. By default (parameter save is set to TRUE) the plot is saved in the working directory.
file.name	Optional. If argument save was set to TRUE specifies the file name for the plot. The default name is "ampliconduo_<date>.<file extension>"
format	Optional. Character specifying the format of the saved file. One of "eps", "ps", "tex", "pdf", "jpeg", "tiff", "png", "bmp", "svg" and "wmf" (windows only). Default format is "jpeg".
...	Optional. Allows to pass other aesthetics.

Author(s)

Anja Lange & Daniel Hoffmann

See Also

[qplot](#), used by [plotAmpliconduo.set](#) to create the plots.

[facet_wrap](#), called for 2-dimensional arrangement of the plots.

[ampliconduo](#), generates the input data.

Examples

```
## loads example data of ampliconduo data frames
data(amplicons)

## plot odds ratio density for amplicon frequencies in ampliconduo data frames
plotORDensity(amplicons)
plotORDensity(amplicons[1:4], binwidth = 0.1, color = "magenta")
plotORDensity(amplicons[[1]], binwidth = 0.1, color = "orange", main = "Sample FU25")
plotORDensity(amplicons[1:4], color = "darkblue", ncol = 2)
```

`site.f`*Sampling Site Short Names*

Description

Character vector with the names of the sampling sites, corresponding to the names used to denote amplicon frequencies in the [ampliconfreqs](#) data.

Usage

```
data(site.f)
```

Format

The format is: chr [1:8] "FU25" "FU28" "FU31.1" "FU31.2" "FU34" "FU37" "UniPond" "BogSoil"

References

Boenigk J, Heider D, Jost S, Lange A, Budeus B, Schilling E, Strittmatter A, Hoffmann D: A high-throughput amplicon sequencing and analysis protocol for comparative analyses of microbial communities (*submitted*)

Examples

```
data(site.f)
data(ampliconfreqs)
ampliconduo(ampliconfreqs[,1:6], sample.names = site.f[1:3])
```

Index

- * **datasets**
 - ampliconfreqs, [7](#)
 - ampliconfreqs.long, [8](#)
 - amplicons, [9](#)
 - site.f, [19](#)
- * **hplot**
 - plotAmpliconduo, [13](#)
 - plotAmpliconduo.set, [15](#)
 - plotORdensity, [17](#)
- * **htest**
 - ampliconduo, [5](#)
 - discordance.delta, [9](#)
- * **list**
 - accepted.amplicons, [4](#)
- * **manip**
 - filter.ampliconduo, [11](#)
 - filter.ampliconduo.set, [12](#)
- * **package**
 - AmpliconDuo-package, [2](#)

accepted.amplicons, [3](#), [4](#), [11](#), [13](#)
AmpliconDuo (AmpliconDuo-package), [2](#)
ampliconduo, [2](#), [4](#), [5](#), [9–18](#)
AmpliconDuo-package, [2](#)
ampliconfreqs, [7](#), [19](#)
ampliconfreqs.long, [7](#), [8](#), [9](#)
amplicons, [9](#)

discordance.delta, [3](#), [6](#), [9](#)

facet_wrap, [16–18](#)
filter.ampliconduo, [3](#), [4](#), [11](#), [13](#)
filter.ampliconduo.set, [3](#), [4](#), [12](#)
fisher.test, [5](#), [6](#)

p.adjust, [5](#), [6](#)
plotAmpliconduo, [2](#), [6](#), [13](#), [17](#)
plotAmpliconduo.set, [2](#), [6](#), [15](#), [15](#), [17](#), [18](#)
plotORdensity, [2](#), [17](#)

qqplot, [15–18](#)

site.f, [7](#), [8](#), [19](#)
xtable, [10](#)